

LIEBHERR





EASY GUIDE INSTRUCTIONS LIEBHERR LABORATORY REFRIGERATORS

1. Switching the appliance on and off

Connect the appliance to the mains - the display reads OFF.

Keep the ON/OFF button  pressed for approx. 5 seconds - the display reads ON.

2. Setting the temperature

- Press  button for 1 second. The temperature display flashes.
- To increase the temperature (warmer): press button .
- To reduce the temperature (colder): press button .
- Press button  again.


The desired temperature setting is saved.

3. Setting the real time clock

The real time clock is preset (CET). Other time zones or summer/winter time must be adjusted manually.

Press  for 5 seconds. Display = *rt1*

 Display = *tc*

 Display = *4 16* Year 2016


 Display = *10* Set the year by pressing the  buttons.

 = save new setting

 Display = *07* Month (1-12)

 Display = *7* Set the month by pressing the  buttons.

 = save new setting

 Display = *01* Day (1-31)


 Display = *1* Set the day by pressing the  buttons.

 = save new setting

 Display = **03** Days of the week (1 = Monday, 7 = Sunday)


 Display = Set the day of the week by pressing the   buttons.

 = save new setting

 Display = **12** Hour (0-23)

 Display = **12** Set the hour by pressing the   buttons.

 = save new setting

 Display = **48** Minute (0-59)

 Display = **48** Set the minutes by pressing the   buttons.

 = save new setting

Press  for 5 seconds.

The electronic control system will automatically switch back to normal operating mode.

4. Calling up stored alarm events and reading the temperature progression

 Display = **HA_n**

Scroll through the list using  or .

HA_n Number of temperature alarms

HA Last temperature alarm

HA 1 Last temperature alarm but one

HA2 Temperature alarm before

HF_n Number of power cuts

HF Last power cut

HF 1 Last power cut but one

HF2 Power cut before **HF 1**

r_t Period in hours in which the maximum and minimum interior temperatures were measured

r_H Maximum (highest) measured temperature

r_L Lowest measured temperature

Select the required item using the  button. Press this button again to return to the list.


Note: You can exit the menu at any time by pressing  for 5 seconds. If no button is pressed within 60 seconds, the electronic control system switches back automatically.

5. Resetting the recorded temperature progression r_t

Complete the following steps if you wish to reset the value saved for r_t in the previous section to 0.

 Display = **HAn**

Press the button  or  until r_t appears in the display.

 Display = **0-999**

Press  for 5 seconds. Display = **rE5**.

The values for r_H and r_L (highest and lowest measured interior temperature) are then reset to the current interior temperature.

Press  for 5 seconds.

The electronic control system will automatically switch back to normal operating mode.